

Kylon® Hooked fabric



Multi-Patented Medical Fabric

Sickle-shaped flexible nylon hooks function as a soft brush when upright, and transform into micro-curettes when pressed firmly against tissue

Minimally invasive, single-use, disposable, sterile multi-purpose “brush-curettes” for tissue sampling and debridement



SFT-1100

SoftBiopsy®+D for debridement and wound base biopsy



SFT-4500

Soft K-Rette® for crevice/web space wound base debridement or tissue sampling



SFT-4100

Soft K-Bride® for minimally invasive and clinically effective wound debridement and tissue removal



Integrated centimeter markings on one handle side for wound depth, width, and length documentation

Coding and Reimbursement - CPT Coding
Low Cost Disposable - Misc. HCPCS Code

Existing wound biopsy and debridement codes are reimbursed with use of Kylon® devices - no change in billing for providers or institutions

Wound Biopsy and Debridement Codes

CPT -1102/3 - Wound Tangential (curettage) Biopsy
CPT-97602- Institutional non-selective debridement coding
CPT 9759 7/8 - Debridement (e.g., high pressure waterjet with/without suction, sharp selective debridement with scissors, scalpel and forceps), open wound (e.g., fibrin, devitalized epidermis and/or dermis, exudate, debris, biofilm), including topical application(s), wound assessment, use of whirlpool, when performed and instruction(s) for ongoing care, per session; total wound(s) surface area; each additional 20 sq cm, or part thereof (List separately in addition to code for primary procedure)

CPT 11042/5 - Debridement, subcutaneous tissue (includes epidermis and dermis, if performed); each additional 20 sq cm, or part thereof (List separately in addition to code for primary procedure)

CPT 11043/6 - Debridement, muscle and/or fascia (includes epidermis, dermis, and subcutaneous tissue, if performed); each additional 20 sq cm, or part thereof (List separately in addition to code for primary procedure)

CPT 15002/3 - When Kylon® device is used for a surgical procedure to excise wound devitalized tissue or eschar to prepare the base for a skin graft

A code: 4649
Surgical Device

Scan for codes



SFT-4200

Kylon® fabric

Soft K-Cot® for gentle fingertip wound debridement and tissue sampling

4095 E. La Palma Avenue, Suite N, Anaheim, CA 92807
Toll Free: 888-235-2275 | histologicswc.com
For ordering information, email support@histologicswc.com